Starla

For Straight Stairs

Stannah

The Stairlift People

The Model 600 with Starla chair is our premium stairlift for straight stairs. With a high, rounded back and gently tilted seat, the Starla is the most comfortable stairlift we've ever made. A button on the armrest folds the footrest up and down, making it easy to get in and out of the seat.

Choose from five upholstery colors, which are also available trimmed with light- or dark wood. The slim-line rail has a subtle bronze sheen, and sits neatly at the edge of the stairway. Quickly installed, the Starla 600 is an exceptional choice which will save you the trouble and expense of moving to a new home.



The Stannah 600 with Starla chair is designed with upholstery choices to blend into your home, and runs along a slim rail attached to your stairs.



Easily folded without bending over, the Starla chair leaves room for walking past. Use remotes on either landing to call or send your stairlift.



With sensible safety features, this stairlift is a practical way to enhance your peace-of-mind.



When swiveled to face the top landing, the Starla chair forms a barrier across the top of the stairs.

Starla 600 Key Features

Comfortable and safe

The generously padded seat has a high back for comfort, and the seat tilts back slightly for added security.

Plain or wood-trimmed upholstery

Choose fabric or vinyl upholstery, and add dark or light wood trim.

Easy-to-use controls -

The slimline control allows you to use any part of your hand to operate the stairlift.

Seat belt

A retractable, seat belt fastens and unfastens with one hand, and is mounted on the armrest, within easy reach.

Seat sensor

Stops the stairlift from moving until you are securely seated.

Safety sensors

Stop the stairlift if it meets an obstruction.

Rated load

One person, 300lbs



Model shown: Starla chair with Beige vinyl upholstery.

Powered footrest

Folds the footrest up or down. This makes it easier to get in and out of the chair, and avoids the need to bend when folding the footrest.

Arm isolation

Prevents the control on the chair from working until the downside-armrest is lowered

On/Off switch

A pulsing blue light indicates that the switch is off.

Swivel levers

Use the levers, or optional one or two-way powered swivel, to get off your stairlift in the safest position at the top or bottom of your stairs.

Easy folding

Armrest, seat and footrest fold neatly away.

Key Lock

A key lock prevents unauthorized use.

Additional Upholstery Colors



Brown (vinyl)



Stone (woven)



Red (woven)



Green (woven)



All colors are available with dark- or light-wood trim

Your local provider:



(855) DME-LIFT (855-363-5438) info@dmelift.com www.dmelift.com

Serving Illinois - Wisconsin - Indiana

Stannah

Stannah Stairlifts Inc.

Head Office & Showroom: Tel: (508) 520-7878 20 Liberty Way, Suite A, Franklin, MA 02038

E-Mail: Stannah@stannah.com Website: www.Stannah-Stairlifts.com